

節亮康復產品有限公司產品手冊

## **Health Product**

Your Smart Health Care Solution Partners
Health care product | Smart home care | Air purification

## **Eating Improver**

**Z5** 

Ziglite Smart Health Care Products Co., Ltd.



## Company Profile 公司簡介

#### ABOUT US 關於我們

Ziglite Smart Health Care Products Co., Ltd. Incorporated in 2018, we work with our partners to bring their products to the local market. Mainly engaged in the management of health care products, the main customers are hospitals and clinics, homes for the elderly and home patients.

We are committed to the development of excellent health care products, dedicated to the needs of a wide range of different customers, get the best care in the family treatment process, and greatly improve the effect of rehabilitation. Get the best care during the family treatment process and greatly improve the healing effect. We will provide quality products with stable quality and stable performance for customers.


We specialize in a wide range of products including: anti-virus and air freshening, incontinence care, mobility aids, decompression supplies, etc.

# what's inside...

公司於2018年成立,我們與合作夥伴一起把他們的產品推向本地市場。以 經營康護用品為主,主要客戶為醫院及診所,安老院及家居病人等。

我們致力於發展卓越的康護醫療產品,竭誠照顧廣大不同客戶的需要,在家庭治療過程中得到最好的護理,大大提高康復治療效果。在家庭治療過程中得到最好的護理,大大提高康復治療效果。務求為客戶提供品質優越及性能穩定的醫療產品。

我們專營的產品種類繁多,包括:抗病毒及空氣清新,失禁護理,助行用品, 減壓用品等。



## **Food Aid Tableware**

## 助食餐具





## GYENNO-Intelligent Anti Shake Spoon

### 睿餐智能防抖勺

#### Model: DC

#### **Features:**

- Portable and replaceable meal head, easy to remove the meal head, cleaning without obstacles
- Small and light, easy to carry
- Effectively assist in shaking eating, eating data, and automatic identification and correction of meal heads
- Medical Tritan shell + medical ultra-thin silica gel, nonirritating, non-toxic, and non-allergic
- 可攜式可替換餐頭,輕鬆拆卸餐頭,清洗無障礙
- 小巧輕便,隨行隨帶
- 有效輔助抖動進食、進食數據、餐頭自動識別矯正
- 醫用Tritan外殼+醫用超薄矽膠,無刺激性、無毒性、無過敏反應



#### **Specifications**

Product Model 產品型號	DC
Product Material 產品材質	Tritan shell/light and soft thin silicone Tritan外殼/輕柔軟薄矽膠
Product function 產品功能	Feeding aid 進食輔助
Interface 介面	USB 2.0
Battery Capacity 電池容量	700mAh
Charging Time 充電時長	2.5H
Standby Time 待機時長	180min
Packing Size 包裝尺寸	17.8*5.5*5.2cm
Package Weight 包裝重量	800g

#### Instructions

- ♦ Long press the top of the spoon handle for 3 seconds to start (red is low battery, please charge)
- ♦ Intelligently correct the direction of the spoon head; tap the function button, the position of the spoon head automatically rotates 360, which is convenient for rolling noodles and digging food
- ◆ 長按勺柄頂部3秒啟動(紅色為低電量,請充電)
- ◆ 智能矯正勺頭朝向,輕點功能鍵,勺頭位置自動旋轉 360度,便於卷麵條、挖取食物

## Bendable Stainless Steel Spoon

## 可彎曲不銹鋼助食勺子

#### Model: JAQ

#### Features:

- The material is stable and healthy and resistant to high temperature, head 304 stainless steel, adjustable angle
- Adjustable hand strap, food grade silicone material, safe and healthy
- Sanitary, detachable, easy to wash and convenient to carry
- Applicable people: hand spasm, finger deformation, finger trembling, loss of grip strength, the elderly, children, people with hand injury or disability
- 材料穩定健康,抵抗高溫,頭部304不鏽鋼,可調節角度
- 可調節手帶,食用級矽膠材質,安全健康
- 衛生、可拆卸,易洗、攜帶方便
- 適用人群:手部曲痙攣、手指變形、手指顫抖、握力喪失者,老人、兒童、手部受傷或殘疾人士



#### **Specifications**

Product Model 產品型號	JAQ
Product Material 產品材質	304 stainless steel, food grade silicone 304不銹鋼,食用級矽膠
Product function 產品功能	Feeding aid 進食輔助
Size 尺寸	24*5*3.3cm
Net Weight 淨重	102g
Packing Size 包裝尺寸	26*5.5*3.5cm
Package Weight 包裝重量	126g

- ♦ Open and take out the appliance components, adjust and use as needed.
- ◆ 打開並將器具組件取出,根據需要調整使用。



## Flexible Stainless Steel Fork

### 可彎曲不銹鋼助食叉子

#### Model: JAR

#### Features:

- The material is stable and healthy and resistant to high temperature, head 304 stainless steel, adjustable angle
- Adjustable hand strap, food grade silicone material, safe and healthy
- Sanitary, detachable, easy to wash and convenient to carry
- Applicable people: hand spasm, finger deformation, finger trembling, loss of grip strength, the elderly, children, people with hand injury or disability
- 材料穩定健康,抵抗高溫,頭部304不鏽鋼,可調節角度
- 可調節手帶,食用級矽膠材質,安全健康
- 衛生、可拆卸,易洗、攜帶方便
- 適用人群:手部曲痙攣、手指變形、手指顫抖、握力喪失者,老人、兒童、手部受傷或殘疾人士



#### **Specifications**

Product Model 產品型號	JAR
Product Material 產品材質	304 stainless steel, food grade silicone 304不銹鋼,食用級矽膠
Product function 產品功能	Feeding aid 進食輔助
Size 尺寸	24.5*5*3.3cm
Net Weight 淨重	103g
Packing Size 包裝尺寸	26*5.5*3.5cm
Package Weight 包裝重量	126g

#### Instructions

- Open and take out the appliance components, adjust and use as needed.
- ◆ 打開並將器具組件取出,根據需要調整使用。

## Anti-shifting And Anti-sprinkling Bowl

### 防移位防灑碗

#### Model: JAP

#### Features:

- The protective design can prevent the flat bowl from scalding or spilling food
- Soft bottom suction cup, the bottom of the bowl with soft silicone
- Press exhaust adsorption, so that the bowl fixed on the table is not easy to slide
- Food grade PP material, safe and assured
- Suitable for: elderly and disabled people who can't carry bowl, patients with hand stiffness, stroke, contracture, children, etc
- 護擋設計,限避免平沿碗使用時造成燙傷或溢出食物
- 軟底吸盤,碗底採用柔軟的矽膠
- 按壓排氣吸附,使碗固定在桌面不易滑動
- 食品級pp材料,安全放心
- 適用於:無法端碗的老年及殘疾人,手僵硬、中風、攣縮病人, 兒童等



#### **Specifications**

Product Model 產品型號	JAP
Product Material 產品材質	Food grade PP plastic + Non-toxic grade silicone 食用級PP塑膠+無毒級矽膠
Product function 產品功能	Feeding aid 進食輔助
Size 尺寸	12.5*8*6cm
Net Weight 淨重	75g
Packing Size 包裝尺寸	18.2*12.5*8.2cm
Package Weight 包裝重量	117g

- ◆ The bottom bracket is seamlessly connected during installation, and it is convenient to disassemble. It can be used separately when it is not required to be fixed. The suction cup can be removed for heating in a microwave oven.
- ◆ 安裝時底托無縫銜接,拆卸方便,無需固定時可分離使用,拆下吸盤可用於微波爐內加熱使用。



## Self-elastic Anti-slip Chopsticks

### 自彈式防滑夾菜筷子

#### Model: JAS

#### Features:

- The latch design can automatically open the chopsticks and close the chopsticks after the food is clamped
- Fit the finger shape to match the palm design
- Curved design, in line with human body force balance
- Applicable people: people with hand cramps, deformed fingers, trembling fingers, the elderly, children, people with hand injuries
- 插銷設計,可以自動打開筷子,夾食後手力閉合
- 採用貼合指形,符合手掌設計,滿足多種人握筷習慣
- 弧形設計,符合人體受力平衡,增進送食的舒適感和不同的角度
- 適用人群:手部曲痙攣、手指變形、手指顫抖者,老人、兒童、手 部受傷或殘疾人士



#### **Specifications**

Product Model 產品型號	JAS
Product Material 產品材質	Plastic + Bamboo 塑膠+竹子
Product function 產品功能	Feeding aid 進食輔助
Size 尺寸	22*5.5*3.3cm
Net Weight 淨重	50g
Packing Size 包裝尺寸	26*5.5*3.5cm
Package Weight 包裝重量	64g

#### Instructions

- ◆ Open and take out the appliance components to use
- ♦ Bamboo chopsticks head can be removed and replaced

06

- ◆ 打開並將器具組件取出,即可使用
- ◆ 竹製筷頭,可去下自行替換

## Stainless Steel Spoon

#### 不鏽鋼助食勺子

#### Model: JAT

#### Features:

- The material is stable and healthy, can be disinfected and resistant to high temperature
- Sanitary, detachable, easy to wash and carry
- The spoon head is 304 stainless steel, which can rotate 360 degrees. It is suitable for different angles and is more comfortable to use
- Adjustable Bracelet, food grade silicone material, safe and healthy
- Applicable crowd: hand spasm, finger deformation, finger shaking, grip loss, elderly, children, hand injury or disabled
- 材料穩定健康,可消毒,抵抗高溫
- 衛生、可拆卸,易洗、攜帶方便
- 勺頭304不鏽鋼,可360度旋轉,適合於不同角度,使用更舒適
- 可調節手帶,食用級矽膠材質,安全健康
- 適用人群:手部曲痙攣、手指變形、手指顫抖、握力喪失者,老人、兒童、手部受傷或殘疾人士





24cm

#### **Specifications**

Product Model 產品型號	JAT
Product Material 產品材質	304 stainless steel, food grade silicone, TPE soft rubber, ABS 304不銹鋼、食用級矽膠、TPE軟膠、ABS
Product function 產品功能	Feeding aid 進食輔助
Size 尺寸	24*5*3.3cm
Net Weight 淨重	84g
Packing Size 包裝尺寸	26*5.5*3.5cm
Package Weight 包裝重量	100g

- ♦ Open and take out the appliance components, adjust and use as needed.
- ◆ 打開並將器具組件取出,根據需要調整使用。



## Stainless Steel Food Fork

#### 不鏽鋼助食叉子

#### Model: JAU

#### Features:

- The material is stable and healthy, can be disinfected and resistant to high temperature
- Sanitary, detachable, easy to wash and carry
- The spoon head is 304 stainless steel, which can rotate 360 degrees, it is suitable for different angles
- Adjustable Bracelet, food grade silicone material, safe and healthy
- Applicable crowd: hand spasm, finger deformation, finger shaking, grip loss, elderly, children, hand injury or disabled
- 材料穩定健康,可消毒,抵抗高溫
- 衛生、可拆卸,易洗、攜帶方便
- 勺頭304不鏽鋼,可360度旋轉,適合於不同角度
- 可調節手帶,食用級矽膠材質,安全健康
- 適用人群:手部曲痙攣、手指變形、手指顫抖、握力喪失者,老人、兒童、手部受傷或殘疾人士





24.5cm

#### **Specifications**

Product Model 產品型號	JAU
Product Material 產品材質	304 stainless steel, food grade silicone, TPE soft rubber, ABS 304不銹鋼、食用級矽膠、TPE軟膠、ABS
Product function 產品功能	Feeding aid 進食輔助
Size 尺寸	24.5*5*3.3cm
Net Weight 淨重	78g
Packing Size 包裝尺寸	26*5.5*3.5cm
Package Weight 包裝重量	94g

#### Instructions

- ♦ Open and take out the appliance components, adjust and use as needed.
- ◆ 打開並將器具組件取出,根據需要調整使用。


## **Chopsticks For The Elderly**

#### 老人康復助食夾菜筷子

#### Model: JAV

#### Features:

- The latch design can automatically open chopsticks and close them by hand after clamping
- Adopt the fitting finger shape, which conforms to the palm design and meets the habit of holding chopsticks by many people
- The arc design conforms to the force balance of the human body and enhances the comfort and different angles of feeding
- Applicable people: people with hand spasm, finger deformation and finger trembling, the elderly, children, people with hand injury or disability
- 插銷設計,可以自動打開筷子,夾食後手力閉合
- 採用貼合指形,符合手掌設計,滿足多種人握筷習慣
- 弧形設計,符合人體受力平衡,增進送食的舒適感和不同的角度
- 適用人群:手部曲痙攣、手指變形、手指顫抖者,老人、兒童、手 部受傷或殘疾人士



#### **Specifications**

JAV
Plastic + Bamboo 塑膠+竹子
Feeding aid 進食輔助
22*5.5*3.5cm
39g
26*5.5*3.5cm
54g

- Open and take out the appliance components to use
- Bamboo chopsticks head can be removed and replaced
- ◆ 打開並將器具組件取出,即可使用
- ◆ 竹製筷頭,可取下自行替換



## Liquid Feeding Straw Aid

#### 流食吸管助食器

#### Model: JEZ

#### Features:

- The material is resistant to falling, easy to remove and replace, and easy to clean after use
- The cup cover can rotate 360 degrees, which is more convenient for nursing
- Clear scale and clear drinking capacity
- Applicable population: inpatients, postoperative patients, pregnant women, long-term bedridden elderly and other people in need of care
- 材質耐摔,拆換便捷,使用後易清洗
- 杯蓋可以360度旋轉,護理更方便
- 刻度清晰,飲用量一目了然
- 適用人群:住院病人、術後病人、產婦、長期臥床老年人等需要 護理人群



#### **Specifications**

Product Model 產品型號	JEZ
Product Material 產品材質	PP
Product function 產品功能	Feeding aid 進食輔助
Capacity 容量	350ml
Pipette Length 吸管長度	28cm
Size 尺寸	8.5*8.5*10cm
Net Weight 淨重	161g

#### Instructions

- ◆ 飲用液體時,可使用軟管接頭杯蓋採用透明材質,方便 觀察飲吸管可隨意彎曲,拆卸方便
- ◆ 飲用流食時,可使用流食接頭杯蓋採用透明材質,方便 觀察防漏設計手抖也可以放心使用
- ◆ 去掉杯蓋可作為日常杯使用大容量且耐磨抗摔,適合兒 童及老人使用

## Feeding Syringe

#### 直咀餵食針筒

#### Model: JFA

#### Features:

- Good sealing without leakage, high contact density between piston head and pipe wall, no gap left
- Impermeable jacket: easy to observe liquid level and bubbles, thickened and durable
- Core rod: anti-skid structure to prevent the core rod from slipping accidentally
- Scope of use: feeding, nasal feeding, washing, perfusion and other multi-purpose
- 密封良好不洩露,活塞頭與管壁接觸密度高,不留空隙
- 透阻外套:便於觀察液面及氣泡,加厚耐用
- 芯杆:防滑脫結構,防止芯杆意外滑落
- 使用範圍:餵食、鼻飼、沖洗、灌注等多用途



#### **Specifications**

Product Model 產品型號	JFA
Product Material 產品材質	PP
Product function 產品功能	Feeding aid 進食輔助
Capacity 容量	20ml,60ml,100ml
Size 尺寸	20ml 14.5*2*4cm 60ml 18.5*5*4cm 100ml 20*7*5cm
Net Weight 淨重	20ml 13g 60ml 32g 100ml 48g

- ◆ 注意事項:
- 1) 請勿使用金屬清潔球等尖銳物品進行清洗
- 2) 使用後請及時洗淨晾乾
- 3) 請勿長期置於高溫場所,以免變形
- 4) 請勿接觸明火及腐蝕性液體



## Feeding Syringe

#### 餵食針筒

#### Model: JFM

#### Features:

- The seal is good without leakage, the contact density between piston head and pipe wall is high, and no gap is left
- Impermeable jacket: easy to observe liquid level and bubbles, thickened and durable
- Core rod: anti-skid structure to prevent the core rod from slipping accidentally
- Scope of use: feeding, nasal feeding, washing, perfusion and other multi-purpose
- 密封良好不洩露,活塞頭與管壁接觸密度高,不留空隙
- 透阻外套:便於觀察液面及氣泡,加厚耐用
- 芯杆:防滑脫結構,防止芯杆意外滑落
- 使用範圍:餵食、鼻飼、沖洗、灌注等多用途



#### **Specifications**

Product Model 產品型號	JFM
Product Material 產品材質	PP
Product function 產品功能	Feeding aid 進食輔助
Capacity 容量	3ml、5ml、10ml
Size 尺寸	3ml 8.5*2.5*1.5cm 5ml 9*3*2cm 10ml 10*3.5*2cm
Net Weight 淨重	3ml 3g 5ml 6g 10ml 7g

#### Instructions

- ◆ 注意事項:
- 1) 請勿使用金屬清潔球等尖銳物品進行清洗
- 2) 使用後請及時洗淨晾乾
- 3)請勿長期置於高溫場所,以免變形
- 4) 請勿接觸明火及腐蝕性液體



#### Ziglite Smart Health Care Products Co., Ltd. 節亮康護產品有限公司

#### Hong Kong Office Addres

Room B3, 18/F Bonsun Industrial Building, 364-366 Sha Tsul Road, Tsuen Wan, Hong Kong

Office Hour: Mon - Fri: 9:30 - 5:30pm

Email: info@ziglite.com Office: +852 3107 7500 Fax: +852 3544 0462 Whatsapp: +852 6827 2010





